

# Laboratory Freezers

# B Medical Systems | F Range

Laboratory Freezers are devices intended for the safe storage samples, specimens, live virus vaccines, cultures, test materials, chemicals, reagents, and other laboratory preparations at temperatures below -20°C.

Compliant to DIN 13277 | Medical Device according to MDR (EU) 2017/745, Class I or 21CFR Part 862.2050, Class I











SAVING LIVES
THROUGH RELIABLE
AND INNOVATIVE
TECHNOLOGY

Medical Refrigeration



# Safety Standards **B Medical Systems**

**The Safety Standards developed by B Medical Systems** define certain significant technical features of a product. These ensure the safe storage of the preparations as well as setting the highest standards of safety for the user.

PRECISION LINE	PREMIUM LINE	
		B Medical Systems Electronics
		B Medical Systems digital display
		B Medical Systems 7" full touchscreen display
		Safety door lock (with 2 keys) and key-operated power switch ON/OFF (with 2 keys)
		Power indicator light and digital temperature indicator (display: 0.1 digits)
		Controlled fan cooling system for constant temperature and even temperature distribution across the entire refrigerating chamber
		Automatic switch-off of the evaporator fans when the door opens
		Self-contained alarm system with integrated battery takes over the alarm function and temperature value measurements in case of system failure for at least 48 hours
		Acoustic and visual alarm signal in case of temperature alarm and system failure
		The alarm history function on the electronic stores all the relevant values during a temperature alarm, such as: min., max., average temperature and alarm duration
		Remote transmission alarm signal (via potential-free contact) in case of temperature alarm (change-over contact)
		Door opening alarm (visual / acoustic)
		Designed and tested for climatic class SN (ambient temperature range +10°C to +32°C)
		Designed and tested for climatic class SN / T (ambient temperature range +10°C to +43°C)
		Interior made from stainless steel
		Additional remote transmission alarm signal (via potential-free contact) in case of system failure (change-over contact)
		Smooth castors for optimum flexibility of movement
	$\circ$	External water cooling
	0	Ambient temperature sensor
		Ethernet interface for the visualization of all operating and control functions (hardware and software settings) via B Connected monitoring software on a peripheral device (computer)
0		B Connected - Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data
0		DCU (Data Communication Unit) - Hardware module monitoring all operating conditions and passing them through to a central database – via local network (in combination with °B Connected)
	0	Integrated remote temperature monitoring device (RTMD), offering real-time worldwide remote monitoring, data access over WEB and GPS position

# 

# Laboratory Freezers

# B Medical Systems | F Range

7 models • Volume 121 > 949 L • Set temperature -41°C / -32°C • Climate class SN | SN/T • Compliant to DIN 13277 | MDR (EU) 2017/745, Class I

In conformity with national and international guidelines, regulations for Medical Devices offering reliability, efficiency and safety at an optimal price.

#### Integrated multifunction monitoring electronics

Integrated B Medical Systems multifunction electronics; with easy access in the door level; offers alarms, central alarm system, set point security, battery backup, full display and many other control features; compatible with "B Connected solution for data monitoring and recording.



#### • Rotomoulded or steel models

 Model F130: One piece material cabinet with a lifetime warranty for zero corrosion, highest quality polypropylene for better insulation and no sharp edges to provide easier maintenance and cleaning

 Models F290-380: Steel cabinets feature lower overall energy consumption, superior temperature distribution and a highly stable air circulation within the cabinet, remarkable sealing between gaskets, long autonomy and holdover times, and an exceptional life-span for the

 Smart ergonomy - Easier and safer handling because of heavier components placed at the bottom



#### High storage capacity

High storage per square meter in terms of net volume, with a very versatile inner volume utility due to multiple user friendly rails allowing modularity.



#### Exclusive integrated electronics handle bar


Laboratory freezers are devices intended

for the safe storage samples, specimens, live virus

vaccines, cultures, test materials, chemicals, reagents, and other laboratory preparations at temperatures below -20°C developed as a result of 40 years of expertise in refrigeration technology.

All functionalities are easily accessible:







#### High quality materials

High quality coating, certified medical devices quality and antibacterial, high quality steel for better longevity and easy hygiene control.



#### Versatile modularity

Excellent storage capacity and modularity - Large choice of shelves and drawers offering modularity for every need.

#### Special features for easy deicing and defrosting

- Evaporator offering better performance than market average, resulting in lesser ice formation with automatic defrost
- Insulated inner doors for significantly lower loss of cooling when open







# **Technical Data**General features















		F13	0		F29	0		F38	0		F40	0		F50	0		F70	0		F90	00	
Gross / Net volume (I)		121/	106		296 /	281		338 /	318		478 / 3	387		634 / 5	514		791/	641		949 /	768	
Set	temperature (preset)					-32°C										-4	1°C					
	temperature (setting range) be adjusted in steps of 0.1°C					-										-41°C	to -20°C					
Pres	set cold / warm alarm limit				-;	37°C / -2	7°C									-46°C	/ -32°C					
Hole	d over time	0.7 h (	-32°C to -2	3°C)	0.5 h (	-32°C to -2	3°C)	0.6 h (	-32°C to -23	3°C)	2.5 h (-	41°C to -18	3°C)	3.0 h (-	41°C to -18	B°C)	3.0 h (	-41°C to -18	3°C)	3.2 h (	-41°C to -1	8°C)
Clin	nate class (ambient temperature range)				SN	(+10°C to	+32°C)								9	SN / T (+1	L0°C to +43	s°C)				
Def	rosting technique					Manua	ıl									Automa	tic (hot ga	ns)				
Refi	rigerant type											R290										
External dimensions H x W x D (mm)		830 x 595 x 695 1750 x 603 x 650		1950 x 603 x 650		1988 x 699 x 1039 1988		1988 x	8 x 845 x 1039 198		1988 >	988 x 992 x 1039		1988 x 1139 x 1039								
Inne	er dimensions H x W x D (mm)	630 x 475 x 470 1525 x 475 x 421		1735 x 475 x 421		1173 x 447 x 738		38	1173 x 740 x 738		1173 x 887 x 738											
Net	weight with standard equipment (kg)	58	60	60	97			106			276			301			319			346		
	Supply voltage (V)	220 -240	220 -240	115- 127	220 -240	220 -240	115- 127	220 -240	220 -240	115- 127	230	220	115 -127	230	220	115 -127	230	220	115 -127	230	220	115 -127
	Frequency (Hz)	50	60	60	50	60	60	50	60	60	50	60	60	50	60	60	50	60	60	50	60	60
v Sai	Power (W)	300	300	360	300	300	360	300	300	360	1300	1300	1300	1300	1300	1300	1300	1300	1300	1300	1300	1300
erg	Energy consumption (kWh/24h)	2.8	2.9	3.1	2.2	1.9	1.8	1.9	1.8	1.9	5.4	5.0	4.9	6.6	5.4	5.2	7.4	5.7	5.3	6.8	6.1	7.6
	Heat emission (Kcal/h)	100	100	111	162	162	150	153	153	153	193	178	176	236	192	186	265	205	190	244	219	272
	Compressor running time (%)	51	51	50	54	54	41	50	50	40	50	40	40	42	40	43	50	39	43	50	43	52
	Noise level (dB(A)) (at 1m height & 1m distance)	45	47	47	36	37	37	37	38	38	47	49	50	44	46	47	45	47	48	45	47	48



## Technical Data Specifications







F130	F290-380	F400-500-700-900
1 130	1 230 300	1 400 300 700 300

	F130	F290-380	F400-500-700-900
nsor	NTC, at 0°C ±	: 0.3°C	PT1000, at 0°C ± 0.1°C
r / 100 ml DOW (Silicon Oil)	NTC, at 0°C ±	: 0.3°C	PT1000, at 0°C ± 0.1°C
unction time of the control ystem failure	12 V - 2.3 AH	I / 48h	12 V - 10.0 AH / 48h
idity at +32°C		≤ 75%	
on (polyurethane)	82 mm PU	65 mm PU	65-82 mm PU
ation (polyurethane)	60-65 mm PU	57-82 mm PU	77-82 mm PU
Inner body / door	Polypropylene (cop	polymer UV)	Stainless steel (1.4301)
Outer body / door	Polypropylene (copolymer UV)	Painted steel	DX51D + Z100 coated
Separate interior doors	-		Anodized aluminium
Drawer	Polycarbonate, tr	ransparent	-
Shelf	-		Wire DIN177, PA11 coated
dical Device Regulation		MDR (EU) 2017/745, Class I (model	s: 220-240 V - 50/60 Hz & 230 V - 50 Hz)
?		2014 / 30 / EU (models: 220-240	V - 50/60 Hz & 230 V - 50 Hz)
	Inction time of the control Instem failure  dity at +32°C  In (polyurethane)  Inner body / door  Outer body / door  Separate interior doors  Drawer  Shelf  dical Device Regulation	nsor  NTC, at 0°C ±  NTC, at 0°C ±	NTC, at 0°C ± 0.3°C  Silicon oil)  Interpretation (polyurethane)  Interpretation

European Medical Device Regulation	MDR (EU) 2017/745, Class I (models: 220-240 V - 50/60 Hz & 230 V - 50 Hz)						
EMC directive	2014 / 30 / EU (models: 220-240 V - 50/60 Hz & 230 V - 50 Hz)						
Low voltage directive	2014 / 35 / EU (models: 220-240 V - 50/60 Hz & 230 V - 50 Hz)						
FDA regulation   Medical Device	21CFR Part 862.2050, Class I (models: 115-127 V - 60 Hz & 220 V - 60 Hz)						
cCSAus certificate	Yes (models: 115-127 V - 60 Hz & 220 V - 60 Hz)						
Energy Star certificate	rgy Star certificate Yes (models: 115-127 V - 60 Hz)						
Clean room classification   ISO 14644-1	Class 6	Class 4					













**Equipment** General

		F130	F290	F380	F400	F500	F700	F900
B Medical	Digital display	•	•	•	-	-	-	-
Systems Electronics	7" full touchscreen display	-	-	-	•	•	•	•
Ergonomic har open/close - heavy	ndle duty door lock mechanism	-	-	-	•	•	•	•
Separate interi in order to minimize	ior doors e the loss of cold air	-	-	-	●2 ○4	<b>2</b>   <b>4</b>	<b>2</b>   <b>4</b>	●2 ○4
Ethernet interf	iace	-	-	-	•	•	•	•
°B Connected -	Monitoring Software	0	0	0	0	0	0	0
DCU - Data Col (in combination wit	mmunication Unit h °B Connected)	$\circ$	$\circ$	$\circ$	-	-	-	-
Integrated rem monitoring dev	note temperature vice (RTMD)	-	-	-	0	0	0	0
Temperature compression to the compression of the c	hart recorder e form of a circular chart recorder)	0	$\circ$	$\circ$	$\circ$	$\circ$	$\circ$	0
Ambient tempe	erature sensor	-	-	-	0	0	0	0
Potential-free of the case of system fa		•	•	•	•	•	•	•
Integrated por (installed by custom	t for external sensor <sub>eer)</sub>	•	•	•	•	•	•	•
External water	cooling	-	-	-	0	0	0	0
Rollers		<b>2</b>	• 2	• 2	-	-	-	-
Smooth castor	s with stabilizers	-	-	-	• 4	• 4	• 4	• 4
Door hinge Rig	ht   Left	• 10	•10	• 10	• 10	•10	•10	•10










2 versions:

Front Cover

with or without






ST-SHELF (STAINLESS STEEL)

### Equipment Internal storage

Standard		F130	F290	F380	F400	F500	F700	F900
Drawer with Fro	ont Cover	<b>2</b>	<b>6</b>	<b>•</b> 7	-	-	-	-
Drawer without	Front Cover	-	<b>1</b>	<b>1</b>				
Bottom Drawe	er	<b>1</b>	<b>1</b>	<b>1</b>	-	-	-	-
Wire Shelf		-	-	-	• 5	<b>5</b>	<b>5</b>	<b>5</b>
Optional								
ST-Drawer with Front Cover		-	-	-	5 (max. per unit)			
ST-Drawer without Front Cover		-	-	-	13 (max. per unit)			
ST-Shelf		-	-	-	13 (max. per unit)			
Accessories								
Dividers Set (H to facilitate the had	l 100 mm) ndling and storage management	2 (max. per unit)	<b>6</b>	• 7	-	-	-	-
S-Rack (H 100 x	W 113 x D 638 mm)	-	-	-	15 (max. per unit)	20 (max. per unit)	25 (max. per unit)	30 (max. per unit)
CT De el	with Wire Shelf	-	-	-	2 (max. per tray)	3 (max. per tray)	3 (max. per tray)	4 (max. per tray)
ST-Rack (H 95 x W 134	with ST-Drawer	-	-	-	2 (max. per tray)	2 (max. per tray)	3 (max. per tray)	4 (max. per tray)
x D 648 mm)	with ST-Shelf	-	-	-	2 (max. per tray)	3 (max. per tray)	4 (max. per tray)	4 (max. per tray)
ST-Cover for Wi	ire Shelf	-	-	-	13 (max. per unit)			





SET

S-RACK



ST-COVER (STAINLESS STEEL) Available for Wire Shelf

ST-RACK (STAINLESS STEEL)

# Full functionalities at a glance

## B Medical Systems | Electronics

**PRECISION LINE** | Multifunctional electronics with digital display and easy access, in the door handle; compatible with °B Connected monitoring solution.



#### • Key power switch • O = OFF / I = ON

- To switch the unit on, turn the key switch to position "I". The green power LED lights up
- After the self-check, the temperature inside the device is displayed




#### Digital display

During normal operation the display shows:

- Left-hand side: current time and internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°
  - Right-hand side: the maximum and minimum inside temperature



#### Function buttons





 ENTER button: use the ENTER-key to confirm the actual selection. This applies to the parameter value. The change of a parameter only becomes effective when confirmed



 MUTE / BACK button: in case of an alarm, the acoustic alarm can be muted for a specified time. In Parameter-mode the same key acts as "BACK-button" to go one step back





#### Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display with the inside measured temperature. If several alarms occur simultaneously, the messages will be shown alternately.



A USB port is available for data export only. Measured values of the display sensor are stored for at least 30 days (ring buffer) together with a time stamp.



#### The B Medical Systems Electronics

offers a wide range of adjustment and diagnostic functionalities as well as additional protection / warning operations (via external alarm operations, histories and individual display signals).

**PREMIUM LINE** | Multifunction electronics with 7" full touchscreen display, integrated at optimal level in the door handle with pre-installed connection allowing exclusive °B Connected monitoring functionalities.



# OFF / ON

- To switch on press the ON/OFF button for 1 second. The green LED will light up
  - Switch off protected by password.
     The door will automatically unlock before powering down



# 0000

#### Touchscreen 7" full display

During normal operation the display shows:

- At the top: current voltage, current time, defrosting, gasket heater and login status (user level)
- Left-hand side: temperature of the regulation sensor (displaying this temperature can be disabled through the settings) and internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°
- Right-hand side: touch buttons to access door unlock, clean, menu, alarm mute and message functions



○ USB port & SD card slot

#### Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display together with the inside measured temperature. If several alarms occur simultaneously, the messages will be shown alternately.



#### Keycard reader

Login / Unlock using an NFC card

#### Touch buttons

- DOOR UNLOCK button: unlocks the door of the device
- CLEAN button: the cleaning button allows the user to disable the touch capability of the screen for a short period of time in order to clean the screen. This button is replaced by the logout button when user is logged in
- MENU button: allows access to the menu screen
- ALARM MUTE button: deactivates the acoustic alarm for a predefined period of time
- ACTIVE ALARMS AND WARNINGS button: during a warning or alarm situation, the color of this button and the corresponding level will change. Information on an active alarm and warning is accessible by pressing the message button













# Optimum control and protection

# B Medical Systems | °B Connected & DCU



#### **\*B CONNECTED | MONITORING SOFTWARE**

Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data.



- Unique monitoring software for the full range of Blood Management Solutions and Medical Refrigeration products
- Web-based interfaces for computers and mobile devices
- Modern design for simple and intuitive use
- Graphic display of temperature curves
- Integrated event and activity history of appliances' components
- Data recording on centralized database for long-term archiving
- Easy setting of specific alarm, via email or SMS alerts
- Generation of reports compiling crucial data and activities
- Temperature and detailed device data export for third-party software
- Important cost advantage compared to a traditional circular chart recorder and its spare parts
- REST API to access raw data directly from the database in read-only mode

#### **KEY BENEFITS:**

- Unique interface for the monitoring of the full range of refrigeration devices
- Centralized database providing data access to entire customer network







# DETAILED DISPLAY

- Optimized design for mobile and touch compatible devices
   Support of multiple
- Support of multiple languages
- Real-time data: Current temperatures, set point, upper and lower limit, door state, compressor runtime
   Log data: Events, alarms within date/time,
- event type, description of the event. Event and alarm confirmation function • Parameter settings of each device: Factory
- settings, real-time clock, language, firmware version, RS485 address, IP address

- > Complete & legally safe documentation of temperature data
- > Comprehensive applications and diagnostic functionalities









#### DCU | DATA COMMUNICATION UNIT (PRECISION LINE)

# Hardware module monitoring all operating conditions and passing them through to a central database – via local network.



- Interface connection of B Medical Systems appliances to an existing network
- Graphical user interface displaying temperature and configurating devices and alarms easily
- Direct connection to Ethernet and serial bus RS485
- Digital in/out (programmable and customer-specific use)
- Recording and storing appliance-relevant data
- Integrated USB port allowing to export archived data
- The DCU combined with the °B Connected software replaces the paper temperature recorder
- All data are recorded and saved in the internal storage of the DCU and backed up in °B Connected database if connected
- Several additional self-sufficient temperature sensors (up to 4 PT1000) may be connected
- Humidity sensor input (4-20 mA)

#### **KEY BENEFITS:**

- Central system collecting relevant temperature data of the appliances and their respective operating conditions
- A number of connection abilities allowing flexible upgrades for individual projects

temperature measurements in one single graph







State-of-the-art technology for the exacting needs of the medical world

Our Global Expertise

Reliable solutions for safe vaccination around the world







#### After Sales support and service

We strive to provide you with the highest standards of service; not only through our selected distributors and partners for all your maintenance and service but also our second line trouble shooting and after sales service. This factory-based group of engineers is there to help our partners and yourself to get the best solution for your cold storage needs.



Safe global blood management:

from collection to transfusion,



#### SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems (formerly Dometic Medical Systems) has more than 40 years' experience in the medical refrigeration sector.

The company, formerly known as Electrolux Medical Systems, was founded in 1979 when the World Health Organization approached Electrolux in Vianden, Luxembourg, to create a solution for the safe storage and transport of vaccines around the world. In 2001, Electrolux Medical Systems became part of the Dometic Group, and was renamed Dometic Medical Systems. Having established a legitimate reputation in the medical equipment industry, the company has also become a global leader in vaccine cold chain.

#### B Medical Systems S.à r.l.

17, op der Hei L - 9809 Hosingen, Luxembourg

Tel.: (+352) 92 07 31-1 Fax: (+352) 92 07 31-300 info@bmedicalsystems.com















Since 2019 B Medical
Systems has been committed
to the UN Global Compact
corporate responsibility
initiative and its principles
in the areas of human rights,
labour, the environment
and anti-corruption.

#### Luxembourg, in the heart of Europe

